

X099 Ligate & Electroporate 2-4Kb + >4Kb cDNA

Purpose: Make sure the larger sized fractions are okay & determine their sizes

Ligation

tube	vector	insert	H ₂ O	45°C 5'	19ul 10x ligase buffer 1ul 100u T4 DNA ligase
1	pcDNA1-NotI/BstXI gel purified X099 1ul = 24ng	—	8ul	✓	1ul
2	"	cDNA 2-4Kb X095 1.5ul = 6ng	6.5ul	✓	1ul
3	"	cDNA >4Kb X095 6ul = 6ng	7ul	✓	1ul

15°C overnight

JUN 6 7 2004

X099 cont'd

CREATION OF EXPRESSION LIBRARY POOLS

OFFICIAL

~~DATA:~~

Prepare bacteria

- 1) — Streak out MC1061/P3 onto LB/Kanamycin (15 µg/ml) plate

Prepare reagents and everything else

- 2) Make LB/Amp (15 µg/ml)/Tet (8 µg/ml) plates (15 cm). Make 1 liter total.
- 3) Sterilize:
 eppendorfs (1/electroporation)
 ddH₂O (500 ml)
 10% glycerol (500 ml)
 pasteur pipets (long ones)
 1 M MgCl₂
 1 M MgSO₄
 1 M glucose
- 4) Put in cold room:
 all the above sterilized stuff except medium stuff
 electroporation cuvettes (0.2 cm gap)
 pipet tips (yellows/blues?)
 also sign up centrifuges

~~DATA:~~

Bacteria

- 8) Start overnight cultures of MC1061/P3 from A. Pearson's plate of MC1061/P3 from Ladin lab
 Pick at least 2 colonies into 3 mls each of LB + Kanamycin (15 µg/ml)
 Also streak each onto half a plate of LB/Amp/Tet (only undesirable revertants should grow)

~~DATA:~~

Get bacteria ready for electroporation

- 9) Put 3 ml of starter culture into 250 ml LB + Kanamycin (15 µg/ml)
-
- Grow until 0.5 - 0.7 O.D.

time	O.D. ₆₀₀
7:25	-
10:00	0.320
10:20	0.497
10:35	0.600 ⇒ on ice

Clean up ligation

- 10) Add TE pH 8.0 to 50 µl
- 11) Add 50 µl Phenol/Chloroform/Isoamly alcohol
Vortex, Spin and recover top aqueous layer
- 12) Add 50 µl TE pH 8.0 to organic layer to backextract
Vortex, Spin and recover top aqueous layer and add to previous aq. layer (total = 100 µl)
- 13) Add
 50 µl 1X LPA ✓
 10 µl 3M NaOAc ✓
 250 µl 100% ethanol
- 14) Put at -80°C 30 min 8:30 *rinse w/ 70% EtOH*
- 15) Spin down at 4°C, remove supe and air dry (don't dry completely) *speed vac* *tube 1 dried completely*
- 16) Just before ready to use, resuspend in 8 µl TE (sterile) *Keep on ice*

Use 2 µl/electroporation
 Freeze rest in DNA box -20°C.

Get bacteria ready for electroporation (everything on ice!)

- 17) Put culture into ice water to chill 15 min (swirl occasionally) 10:35 - 10:50
- 18) Spin down in 1 disposable conical tube, 4°C, 15 min, 4000 rpm (2600 xg), 10:55 - 11:10
- 19) Decant most but not all liquid (leave equal volume liquid as in pellet). Add 5 ml sterile water and resuspend gently with pipet.
- 20) Add 250 ml ice-cold ddH₂O (sterile), spin 15 min, 4°C, 4000 rpm
- 21) Repeat steps 19-20 but spin 20 min.
- 22) Pour off as much supe as possible (you'll lose some bugs), add 10% glycerol to 12 mls, gently resuspend cells and spin 8,000 rpm 30 min 4°C in SS-34 (in Falcon 2059 tube)
- 23) Pour off supe getting rid of almost all liquid (you'll lose some cells). You want it thick. Resuspend in 100 µl 10% glycerol (you want it thick). Used 200 µl had 100 µl left over.

During spin periods set up for electroporation

- 24) Make SOC from SOB
- 25) Put electroporator chamber on ice
- 26) Connect pulse controller to gene pulser (connect in front the red to red and black to black). Next time resuspend in 50 µl? or don't add any liquid at all? yes should get only 3 electros/250 mls
- 27) Set to:
200 ohms
25 µF
2.5 kV
- 28) Get everything else ready (Falcon 2059 with 1 ml SOC each, pasteurs, tips, etc)

Electroporation

- Always do controls: water only (neg. control) and uncut vector (positive control)
- 30) Swirl bacteria with sterile yellow tip. Pipet up 40 µl bacteria to tube #1 on ice. Pipet up and down avoiding generation of bubbles. Let sit 30 sec on ice.
 - 31) With fresh tip take up bacteria and put into cuvette as close to bottom as possible without creating bubbles. Quickly shake hard down to bottom (v. important).
 - 32) Take off cap, put in electroporator chamber, pulse
 - 33) Quickly remove cuvette and add 1 ml SOC. Resuspend with pasteur pipet and transfer to 15 ml round bottom and incubate shaking at 37°C, 90 min. 1pm - 2pm
 - 34) Repeat steps - for each electroporation.
 - 35) Put LB/Amp/Tet plates into hood to dry.
 - 36) Plate out 50 - 100 µl/plate to test for electroporation efficiency.
Use 1:100 of positive control
1:5 of ligation mix
undiluted neg. control
Grow overnight 37°C. Store electroporated bacterial cultures at 4°C up to one week.
 - 37) Count colonies.

electro #	DNA	tau	amt plated	colonies	effic. (col/µg)	for 5000 col per plate vol of undi
1	Lig #1 - 2 µl of 8 µl	4.5	10 µl + 90 µl LB	70	1.2×10^6	-
2	0 #2	4.5		2072	3.5×10^7	24 µl
3	#3	4.5		1300	2.2×10^7	38 µl
4	- (2 µl TE)	4.5	100 µl	13	-	-
5	pcDNA1 + 2.3 Kb insert (1 µl = 0.5 µg) 2 µl	4.5	1 µl + 99 µl LB	553	5.6×10^8	-

electro

- #2 pcDNA1 + 2.4 Kb 24 µl/plate \Rightarrow 42 plates = 210,000 clones assume 75% loss - 48 hrs later 30 µl/plate \Rightarrow 32 plates
- #3 pcDNA1 + 2.4 Kb 38 µl/plate \Rightarrow 26 plates = 130,000 clones 48 µl/plate \Rightarrow 20 plates

X099 cont'd


PCR of colonies to check cDNA sizes

19 μ l 10X PCR buffer
 47.5 μ l T7 primer
 47.5 μ l SP6 primer
 0.95 μ l Tag polymerase
 1.52 mix of dNTPs (each 75mM)
 73.5 μ l H₂O
 190 μ l total

3 μ l dATP
 3 μ l dCTP
 3 μ l dTTP
 3 μ l dGTP

12 μ l dNTPs (freeze in PCR box -20°C)

Aliquot 15 μ l PCR oil/tube
 Aliquot 10 μ l above stock soln/tube
 Flame straight needle, poke colony, then into PCR tube
 PCR 94°C 30 sec \rightarrow 50°C 30 sec \rightarrow 72°C 2 min \times 35 cycle
 Add 1 μ l 10X blue juice
 Load 6 μ l/lane onto 0.9% seakem GT6 agarose minigel



colony	DNA size	miniprep
1	-	2.5
2	0.7	
3	0.8	
4	0.7	
5	2.0	
6	1.0	
7	1.9	
8	-	2.1
9	-	2.8

avg size = 1.6
 med size = 1.9

10	1.2	
11	2.4	
12	-	0.8
13	0.9	
14	0.8	
15	1.0	
16	-	3.0
17	-	
18	0.9	

avg. size = 1.4
 med size = 0.95

X099 cont'd

Mini-pups of clones that did not PCR

Method: Maniatis

- Changes: ① spun twice to get rid of white particulate matter after adding SalI III
- ② Phenol/sevage extracted
- ③ Resuspended in 25 μl TE pH 8


original
cDNA
fraction Colonies2-4Kb {
8
9>4Kb {
12
16
17

Digest

plasmid

~~2 μl~~ NEB2
10X buffer7 μl 10X
Ac-BSA7 μl 10 mg/ml
~~200 $\mu\text{g}/\text{ml}$~~
Erase0.1 μl H_2O + 25 μl Not I + Hind III
50:50 mix+ 7 μl 2 μl each + 8 μl above rxn mix

37°C 2 hrs (3:45-5:45)



colony	size
#1	1.1 + 1.4 = 2.5
8	2.1
9	1.0 + 1.8 = 2.8
12	0.8
16	3.0
17	-

PX503 (1 μl each) -

Not sure if state of DNA is
very different or loading
made a difference in the
way it ran

Plasmid midiprep for cDNA library


preps:

B46 - B53

Day 1


1. Scrape 150 mm plate with 5 mls LB. Transfer to Falcon 2059 15 ml tube on ice.
2. Add another 3 mls LB to plate and scrape again.
3. Take 400 μ l, put into freezer vial, add 100 μ l glycerol and freeze at -150C.
4. Spin rest in SS-34, 9000 rpm, 2 min 4C.
5. Dry pellet as much as possible.
6. Resuspend pellet in 500 μ l ice cold solution I by vigorous vortexing.
7. Add 1 ml fresh solution II (0.2 N NaOH, 1% SDS)
for 100 ml: 1 ml 2 N NaOH
0.5 ml 20% SDS
8.5 mls ddH₂O
Swirl gently until clear. Do not vortex. Leave on ice 10 min.
8. Add 750 μ l solution III (ice-cold). Close tube and mix contents by shaking vigorously several times. Store on ice 5'. A flocculent white precipitate should form.
9. Centrifuge 15', 4°C, 9000 rpm.
10. Recover supe and add equal volume of phenol:chloroform. Mix by vortexing.
11. Spin 9000 rpm, 5'.
12. Add 2 volumes, ethanol r.t., vortex, let stand 5'.
Spin in SS-34 rotor for 15' 9,000 rpm.
13. Carefully aspirate off all liquid (want it dry) and rinse once in cold 70% EtOH. Aspirate again getting any extraneous drops. Air dry 10' minutes.
14. Redissolve in 100 μ l TE plus DNase-free RNase (20 μ g/ml). Vortex briefly. Incubate 37°C, 30'.
Transfer to sterile eppendorf.
15. Quantitate by dilutions onto EtBr plate.

Store at 4°C O.N.
Incubate 37°C 2 1/2 hrs
Freeze



Note:

photograph is
mirror image of
schematic above



Slight variation
in props. I think
5 μ transfection
will be good for
all.

X112

Screen cDNA pools B63-B76, redo B47

OFFICIAL

DEAE dextran transfections of COS M6 cells

materials:

- | | |
|---|---|
| 1. 35 mm dishes. | 5. CMF PBS |
| 2. DMEM with 10% FBS. | 6. DEAE-dextran (10 mg/ml in CMF PBS (autoclaved)). |
| 3. Chloroquine (40 mM in CMF PBS, sterile filtered) | 7. DMSO |
| 4. DNA | 8. cPBS |
| | 9. sterile tips |

method:

day 0 (set up cells).

Set COS M6 cells in 35 mm dishes at 300,000 cells/dish in 2 ml DMEM with 10% FBS

day 1 (transfect):

1. In sterile eppendorf's prepare for each dish add (in order):
 - a) DNA - 500 ng/dish
 - b) add CMF PBS to 190 μ l, vortex
 - c) 10 μ l of 10 mg/ml DEAE-dextran, vortex well (7 seconds)

tube #	plates	DNA	CMF PBS	10 mg/ml DEAE-dextran	positive	maybe
1		cDNA pool B63	10 μ l	180 μ l	1	2
2		B64			4	7
3		B65			1	0
4		B66			1	0
5		B67			2	0
6		B68			0	1
7		B69			0	4
8		B70			2	3
9		B71			2	4
10		B72			0	1
11		B73			4	0
12		B74			12	4
13		B75			2	2
14		B76			2	6
15		B47 - redo from X11			9	4
16		B47 - redo from X11			12	2
17		pcDNA1	0.39		0	1
18		1:5000	9.2 μ l		396	8

2. Rinse cells with 2 mls CMF PBS (37C). Aspirate PBS.
3. Add transfection cocktail from step 1 and distribute evenly. Incubate 30 minutes, distributing liquid every 10 min. 9:50 - 10:20
4. Add 2 ml DMEM 10% FBS + 80 μ M chloroquine and incubate 37C 2.5 hrs. 10:20 - 12:50
5. Aspirate off medium and replace with 1 ml 10% DMSO in DMEM 10% FBS for 2.5 min.
6. Aspirate off and wash once with 2 mls cPBS.
7. Refeed with 2 ml warm DMEM 10% FBS/dish. Incubate overnight.

Day 2

Refeed cells with 2 ml DMEM 10% FBS + 1 mM Nabutyrate pH 7.3.

Day 3

Refeed cells with 0.75 ml DMEM 10% FBS + 1 mM Nabutyrate pH 7.3 + 3 μ g/ml DiI-AcLDL for 5 hrs.14.25 ml DMEM 10% + 14.25 μ l 1M Nabut + 158 μ l DiI-AcLDL prep #48 (0.27 μ g/ml)

plate	DNA	positives	maybes
(1)	B47 pool	7 (10)	5 (6)
2	B48	3	4
3	B49	1	0
4	B50	2	6
5	B51	0	1
6	B52	0	0
7	B53	2	2
8	B54	2	5
9	B55	0	0
10	B56	2	4
11	B57	2	5
12	B58	3	8
13	B59	1	4
14	B60	3	2
(15)	B61	6 (4)	2 (3)
16	B62	1	2
17	pcDNA1	1	2
18	1:5000	430	-

Positives are scored if cells are punctate
parentheses are recounts

X120

Create subpools of 15 colonies of B47.1.8

Purpose: Reduce pool size to approx 15 colonies to narrow the search for the MAC26-1 receptor.

Transform competent MC1061/p3 (Q.G. purple dot) usual procedure but didn't incubate on ice 30'; just heat-shocked 37° 5' right away. Still worked

Plated 5ul - count 930 (B47.1.8)
4 (no DNA)

Took 3.2ul transformed bugs + 1.9 mls LB
plated 50ul/plate

Counted plates

	<u>colonies</u>
B47.1.8.1 -	19
.2 -	15
.3 -	19
.4 -	23
.5 -	15
.6 -	13
.7 -	17
.8 -	21
.9 -	15
.10 -	33
.11 -	18
.12 -	17
.13 -	22
.14 -	26
.15 -	11
.16 -	16
.17 -	17
.18 -	15

	<u>colonies</u>
.19 -	25
.20 -	14
.21 -	19
.22 -	16
.23 -	23
.24 -	23
.25 -	25
.26 -	16
.27 -	19
.28 -	23
.29 -	26
.30 -	22
.31 -	19
.32 -	18
.33 -	17
.34 -	20
.35 -	19
.36 -	31

1st 24 pools
452
plate 18.8

all 36 pools
717
19.9

IXIZO cont'd

Plasmid midiprep for cDNA library

preps:

B47.1.8.1 - B47.1.8.24

Day 1

1. Scrape 100 mm plate with 2 mls LB. Transfer to Falcon 2059 15 ml tube on ice.
2. Add another 2 mls LB to plate and scrape again.
3. Spin in SS-34, 9000 rpm, 2 min 4C.
5. Dry pellet as much as possible.
6. Resuspend pellet in 300 μ l ice cold solution I by vigorous vortexing.
7. Add 0.6 ml fresh solution II (0.2 N NaOH, 1% SDS)
for 100 ml: 1 ml 2 N NaOH
0.5 ml 20% SDS
8.5 mls ddH₂O
Swirl gently until clear. Do not vortex. Leave on ice 10 min.
8. Add 450 μ l solution III (ice-cold). Close tube and mix contents by shaking vigorously several times. Store on ice 5'. A flocculent white precipitate should form.
9. Centrifuge 15', 4°C, 9000 rpm.
10. Recover supe and add equal volume (1.2 ml) of phenol:chloroform. Mix by vortexing.
11. Spin 9000 rpm, 5'.
12. Add 2 volumes (2.5) ethanol r.t., vortex, let stand 5'. overnight (1-16) or 6 hrs (17-24)
Spin in SS-34 rotor for 15' 9,000 rpm.
13. Carefully aspirate off all liquid (want it dry) and rinse once in cold 70% EtOH. Aspirate again getting any extraneous drops. Airdry 10' minutes.
14. Redissolve in 50 μ l TE plus DNase-free RNase (20 μ g/ml). Vortex briefly. Incubate 37C, 2 hr.
Transfer to sterile eppendorf.
15. Quantitate by dilutions onto EtBr plate.

X121 Screen subpools B47.1.8.1-B47.1.8.24;
compare CD36 with B47.1.8

DEAE dextran transfections of COS M6 cells

Method:

Day 0 (set up cells)

Set COS M6 cells in 35 mm dishes at 300,000 cells/dish in 2 ml DMEM with 10% FBS

Day 1 (transfect)

1. In sterile eppendorfs prepare for each dish add (in order):

- DNA - 500 ng/dish
- add CMF PBS to 190 μ l, vortex
- 10 μ l of 10 mg/ml DEAE-dextran, vortex well (7 seconds)

tube #	plates	CMF PBS	10 mg/ml DEAE-dextran	Results
1	B47.1.8.1	10 μ l	180 μ l	10 μ l
2	2			
3	3			
4	4			
5	5			
6	6			
7	7			
8	8			
9	9			
10	10			
11	11			
12	12			
13	13			
14	14			
15	15			
16	16			
17	17			
18	18			
19	19			
20	20			
21	21			
22	22			
23	23			
24	24			
25	bsRL (PCR3)	1.5 μ l	570 μ l	30 μ l
26	↓	X	X	X
27	poly I	X	X	X
28	mBSA	X	X	X
29	CD36	1.5 μ l	525 μ l	30 μ l
30	↓	X	X	X
31	poly I	X	X	X
32	mBSA	X	X	X
33	B47.1.8	30 μ l	540 μ l	30 μ l
34	↓	X	X	X
35	poly I	X	X	X
36	mBSA	X	X	X
37	pcDNA1	0.3 μ l	190 μ l	10 μ l
38				

Pool #4 was brightest & w/ more positive cells.

CD36 binds acetylated LDL & is not inhibited by poly I
but is inhibited by small amounts of mBSA
 \therefore CD36 has same properties as MAC 26-1 receptor!

X121 cont'd

2. Rinse cells with 2 mls CMF PBS (37C). Aspirate PBS.
3. Add transfection cocktail from step 1 and distribute evenly. Incubate 30 minutes, distributing liquid every 10 min.
4. Add 2 ml DMEM 10% FBS + 80 μ M chloroquine and incubate 37C 2.5 hrs. 11:50 - 1:20 pm
5. Aspirate off medium and replace with 1 ml 10% DMSO in DMEM 10% FBS for 2.5 min.
6. Aspirate off and wash once with 2 mls cPBS.
7. Refeed with 2 ml warm DMEM 10% FBS/dish. Incubate overnight.

Day 2

Refeed cells with 2 ml DMEM 10% FBS + 1 mM Nabutyrate pH 7.3.

Day 3

Refeed cells with 0.75 ml DMEM 10% FBS + 1 mM Nabutyrate pH 7.3 + 3 μ g/ml DiI-AcLDL for 5 hrs.1.27 mls med + 27 μ l 1M Nabut + 300 μ l DiI-AcLDL (#48 0.27 mg/ml)

2.25 mls + 225 μ l poly I (4 mg/ml) = 400 μ g/ml

2.25 mls + 1.35 m-BSA (3.34 mg/ml) = 2 μ g/ml

9:20 - 2:20 pm

pos cells
4 > 24 > 19

brighter
4 > 24 > 19

1) X122 Create subpools of 1 colony of B47.1.8 and screen

Purpose:

Transform competent (cell) MC1061/P3 (0.6µg purple dot)

1. thaw aliquots (2) of bugs on ice
2. Add 2µl of DNA (or TE) to aliquot
DNA = pool B47.1.8.4
TE = neg. control
3. Carefully pipet up & down twice to mix
4. Heat shock 37°C 5min.
5. Add 200µl LB medium. Shake 1hr 37°C
6. Plate 5µl on 150mm LB Amp/ret plate.

Results: ^{transformation} ~~transformation~~ worked well. Circled 49 apparently single colonies to be picked. Picked each (done by me) *Ana Maria Vranceanu* into 3 ml LB A/T. Grew overnight.

Also had *Ana Maria Vranceanu* repick off same plate & streak onto new plate. (4 indiv/150mm plate).

X122

 cont'd

Plasmid miniprep for cDNA library

preps:

Matrix. 7x7 rows A-F, columns 1-7

14 mini-preps

Day 1

1. Take 200 μ l culture from each tube in a row or column of 7 and put into eppendorf. Store remainder at 4°C.
2. Spin at 12,000 x g for 30 sec in microfuge.
3. Remove medium by aspiration, leaving bacterial pellet as dry as possible.
4. Resuspend pellet in 100 μ l ice-cold solution I by vigorous vortexing.
5. Add 200 μ l fresh solution II (0.2 N NaOH, 1% SDS)
for 2 ml: 0.2 ml 2 N NaOH
0.1 ml 20% SDS
1.7 mls ddH₂O
6. Swirl gently until clear. Do not vortex. Leave on ice 10 min.
Add 150 μ l solution III (ice-cold). Close tube and vortex gently inverted for 5 sec. Store on ice 5'. A flocculent white precipitate should form.
7. Centrifuge 5', 4°C, max speed in microfuge: ^{400 μ l}
8. Recover supe and add equal volume of phenol:chloroform. Mix by vortexing.
9. Spin 2' in microfuge. ^{800 μ l}
10. Add 2 volumes, ethanol r.t., vortex, let stand 2' at r.t.
Spin 5', 4°C max speed in microfuge.
11. Carefully aspirate off all liquid (want it dry) and rinse once in cold 70% EtOH. Aspirate again getting any extraneous drops. Airdry 10' minutes.
12. Redissolve in 10 μ l TE plus DNase-free RNase (20 μ g/ml). Vortex briefly. Incubate 37C, 0.5 hr.
~~Transfer to sterile eppendorf.~~

X122 cont'd Screen matrix

DEAE dextran transfections of COS M6 cells

materials:

1. 35 mm dishes
2. DMEM with 10% FBS
3. Chloroquine (40 mM in CMF PBS, sterile filtered)
4. DNA
5. CMF PBS
6. DEAE-dextran (10 mg/ml in CMF PBS (autoclaved))
7. DMSO
8. cPBS
9. sterile tips

method:

day 0 (set up cells) *6-well*
Set COS M6 cells in *35* mm dishes at 300,000 cells/dish in 2 ml DMEM with 10% FBS

day 1 (transfect)

1. In sterile eppendorfs prepare for each dish add (*isopropyl*):
 - a) DNA - 500 ng/dish
 - b) add CMF PBS to 190 μ l, vortex
 - c) 10 μ l of 10 mg/ml DEAE-dextran, vortex well (7 seconds)

c) 10 μ l of 10 mg/ml DEAE-dextran, vortex well (7 seconds)						Results
tube #	plates	DNA		CMF PBS	10 mg/ml DEAE-dextran	
1	row A	plasmid prep from matrix	5 μ l	18.5 μ l	10 μ l	+ 3 rd brightest
2	B					-
3	C					-
4	D					+ v. faint flou.
5	E					-
6	F					+ 2 nd brightest
7	G					+ brightest row
8	column 1					+
9	2					-
10	3					+
11	4					+ weak
12	5					-
13	6					+ prob.
14	7					+ prob. brighter
15	p.cDNA		0.5	190 μ l		-
16	1:5000 (b5824)		9.3 μ l	180 μ l		-
17						
18						

2. Rinse cells with 2 mls CMF PBS (37C). Aspirate PBS.
3. Add transfection cocktail from step 1 and distribute evenly. Incubate 30 minutes, distributing liquid every 10 min.
4. Add 2 ml DMEM 10% FBS + 80 μ M chloroquine and incubate 37C 2.5 hrs.
5. Aspirate off medium and replace with 1 ml 10% DMSO in DMEM 10% FBS for 2.5 min.
6. Aspirate off and wash once with 2 mls cPBS.
7. Refeed with 2 ml warm DMEM 10% FBS/dish. Incubate overnight.

Day 2

Refeed cells with 2 ml DMEM 10% FBS + 1 mM Nabutyrate pH 7.3.

Day 3

Refeed cells with 0.75 ml DMEM 10% FBS + 1 mM Nabutyrate pH 7.3 + 3 μ g/ml DiI-AcLDL for 5 hrs.

13.5 ml

+ 13.5 μ l 1M Nabut + 150 μ l DiI-AcLDL #48

(0.27 mg/ml)

9:20 - 2:20

Results: Lots! of negatives. Para (a control) test so

X122 (cont'd)

Matrix

	1	2	3	4	5	6	7
A							
B	8	9	10	11	12	13	14
C	15	16	17	18	19	20	21
D	22	23	24	25	26	27	28
E	29	30	31	32	33	34	35
F	36	37	38	39	40	41	42
G	✓ 43	44	45	46	47	48	✓ 49

As stated before,
Plated all 49 out on LB Amp/ret plates
(before results of this exp't were known)

Picked a single colony (hopefully) from
the plates of #7, 28, 42, 43, 45, 46, 48, 49

#42 did not grow, but rest did

Decided to keep things simple so
mini-preped only #43, 45, 48, 49 (see X124)

X124 control

DEAE dextran transfections of COS M6 cells

materials:

1. 35 mm dishes
2. DMEM with 10% FBS
3. Chloroquine (40 mM in CMF PBS, sterile filtered)
4. DNA
5. CMF PBS
6. DEAE-dextran (10 mg/ml in CMF PBS (autoclaved))
7. DMSO
8. cPBS
9. sterile tips

method:

day 0 (set up cells)

Set COS M6 cells in 35 mm dishes at 300,000 cells/dish in 2 ml DMEM with 10% FBS

day 1 (transfect) - Note: cells were way too heavy (used Coulter counter) Its better to

- a) DNA - 500 ng/dish
 - b) add CMF PBS to 190 μ l, vortex
 - c) 10 μ l of 10 mg/ml DEAE-dextran, vortex well (7 seconds)
- use 1 confluent T75 for 18 wells
than use the counter which isn't working.

tube #	plates	CMF PBS	10 mg/ml DEAE-dextran	results
1	847.1.8.4.43	185 μ l	10 μ l	+
2	45			-
3	48			-
4	49			+
5	pcDNA1	190		+
6	PX5R3			+
7				
8				
9				
10				
11				
12				
13				
14				
15				
16				
17				
18				

2. Rinse cells with 2 mls CMF PBS (37C). Aspirate PBS.
3. Add transfection cocktail from step 1 and distribute evenly. Incubate 30 minutes, distributing liquid every 10 min.
4. Add 2 ml DMEM 10% FBS + 80 μ M chloroquine and incubate 37C 2.5 hrs.
5. Aspirate off medium and replace with 1 ml 10% DMSO in DMEM 10% FBS for 2.5 min.
6. Aspirate off and wash once with 2 mls cPBS.
7. Refeed with 2 ml warm DMEM 10% FBS/dish. Incubate overnight.

day 2

Refeed cells with 2 ml DMEM 10% FBS + 1 mM Nabutyrate pH 7.3.

Day 3

Refeed cells with 0.75 ml DMEM 10% FBS + 1 mM Nabutyrate pH 7.3 + 3 μ g/ml DiI-AcLDL for 5 hrs.

6 wells

6.5 mls

1M

0.27 mg/ml

4.85 ml 1M EFT

100

0.45

1M (stock)

B47.1.8.4.49 was chosen as cloned
MAC26-1 receptor &
was renamed phaSR^{III}
(Note: Ana Maria Vincenian
did several exp'ts to show
that the plasmid was a
single one & repeatedly
gave the expected activity)

X133 Sequencing of 5' end of phaSR^{III}
w/ Sequinase Kit

Purpose: Determine if phaSR^{III} encodes
a known protein (e.g. CD36 or
Limp-II)

Followed Kit instructions. Used primers:

T7 - on plasmid pcDNA1

oSA1.3 - in cDNA sequence

5' CTG TCG CTG TCC CCC TTC AG 3'

T7

short read

GGTACCGAGCTCGATCCACTAGTAACGGCCGCCAGTGTGCTCTAAGGCCACCTGCAGGGCTACTG
CTGCTCCGGCCACTGCCTGAGACTCACCTTGCTGGAACGTGAGCCTCGGCTTCTGTCTCTCTG

long read

ACTCACCTTGCTGGAACGTGAGCCTCGGCTTCTGTCTCTCTGTGGCCTCTGTGCTTCTGTGCGT
GTCCCCCTTNAGTCCCTGAGCCCCGCGAGCCCGGGCCGCACACGCGACATGGGCGGCCANNCCAGGG

mix

GGTACCGAGCTCGATCCACTAGTAACGGCCGCCAGTGTGCTCTAAGGCCACCTGCAGGGCTACTG
CTGCTCCGGCCACTGCCTGAGACTCACCTTGCTGGAACGTGAGCCTCGGCTTCTGTCTCTCTGTG
GCCCTGTGCTGCTTCTGTGCTGCTGCTCCCCCTTNAGTCCCTGAGCCCCGCGAGCCCGGGCCGCACACG
GCACATGGGCGGCCANNCCAGGG

OSA3.1

short read

CCGCACACGCGACATGGGCGGCACGGCCAGGGCGCTGGGTGGCGGTGGGGCTGGGCGTCTGTTGGGG
TGCTGTGCGCTGTGCTCGGTGTGGTTATGATCCTCGTGATGCCCTCGCTCATCAAACAGCAGGTAC
TGAAGAATGTC (C?) GCATAGACCCAGCAGCCTGTCTTTGCAATGTGGAAGGAGA

long read

TGTGCTCGGTGTGGTTATGATCCTCGTGATGCCCTCGCTCATCAAACAGCAGGTACTGAAGAATGT
CCGCATAGACCCAGCAGCCTGTCTTTGCAATGTGGAAGGAGATCCCTGTACCCTTCTACTTGTG
CGTCTACTTCTTCGAGGTGGTCAATCCAGCGAGATCCTAAAGGGTGAGAA

mix T7 and OSA3.1

GGTACCGAGCTCGATCCACTAGTAACGGCCGCCAGTGTGCTCTAAGGCCACCTGCAGGGCTACTG
TGCAGGGCTACTGCTGCTCCGGCCACTGCCTGAGACTCACCTTGCTGGAACGTG
AGCCTCGGCTTCTGTCTCTCTGTGGCCTCTGTGCTTCTGTGCTGCTCCCCCTTN
AGTCCCTGAGCCCCGCGAGCCCGGGCCGCACACGCGACATGGGCGGCACGGCC
AGGGCGCTGGGTGGCGGTGGGGCTGGGCGTCTGGGGCTGCTGTGCGCTGTGCT
CGGTGTGGTTATGATCCTCGTGATGCCCTCGCTCATCAAACAGCAGGTACTGA
AGAATGTCCGCATAGACCCAGCAGCCTGTCTTTGCAATGTGGAAGGAGATC
CCTGTACCCTTCTACTTGTCCGTCTACTTCTTCGAGGTGGTCAATCCAGCGAG
ATCCTAAAGGGTGAGAA

BLASTX 1.3.9MP

[Build

Reference: Gish, Warren and David J. States (1993). Identification of protein coding regions by database similarity search. Nature Genetics 3:266-72.
 Altschul, Stephen F., Warren Gish, Webb Miller, Eugene W. Myers, and David J. Lipman (1990). Basic local alignment search tool. J. Mol. Biol. 215:403-410.

Notice: statistical significance is estimated under the assumption that the equivalent of one reading frame in the query sequence codes for protein and that significant alignments will involve only coding reading frames.

Query= TITLE phasr3.seq
 (447 letters)

Translating both strands of query sequence in all 6 reading frames

Database: Non-redundant PDB+SwissProt+PIR+SPupdate+GenPept+GPupdate,
 EDT

96,634 sequences; 27,090,059 total letters.

Searching.....done

Sequences producing High-scoring Segment Pairs:	Reading Frame	High Score	Smallest Poisson Probability P(N)	N	
sp P27615 LIM2_RAT	LYSOSOME MEMBRANE PROTEIN II (L...	+2	114	1.1e-08	1
pir JQ1523 JQ1523	lysosomal membrane 85K sialogly...	+2	109	6.3e-08	1
sp P10284 HM26_MOUSE	HOMEODOMAIN PROTEIN HOX-2.6. >pir ...	-2	61	2.4e-06	2
sp P16671 CD36_HUMAN	PLATELET GLYCOPROTEIN IV (GPV)...	+2	94	1.1e-05	1
gp L06850 HUMCD36B_1	antigen CD36 [Homo sapiens]	+2	94	1.1e-05	1
gp L19658 RATFAT_1	FAT gene product [Rattus norveg...	+2	92	2.3e-05	1
pir A43932 A43932	mucin - human (fragment) 0.0 ...	-1	60	3.8e-05	2
pir B60492 B60492	homeotic protein Hox B4 - human...	-2	57	4.0e-05	2
sp Q01200 PRIA_LENED	PRIA PROTEIN. >pir S23106 S2310...	-1	62	5.6e-05	2
pir S12968 S12968	Acrosin, sperm - Pig #EC-number...	-2	59	6.7e-05	2
gp L23108 MUSCDANTI_1	CD36 antigen [Mus musculus]	+2	88	9.0e-05	1
pir A45106 A45106	mucin - human (fragment) 0.0 ...	-1	60	9.2e-05	2
pir S31976 S31976	Cvx peptide - Rat 0.0 0.0 0.0...	-3	57	0.00012	2
gp Z16406 MOX2A_1	Mox-2 [Mus musculus]	-3	57	0.00012	2
gp Z17223 RNGAKMR_1	Gax peptide [Rattus norvegicus]	-3	57	0.00012	2
sp P13983 EXTN_TOBAC	EXTENSIN PRECURSOR (CELL WALL H...	-2	56	0.00024	2
pir G60110 G60110	repetitive protein antigen 69/7...	-2	81	0.00035	1
gp M14721 MUSFGNAA_1	Mouse epidermal profilaggrin mR...	-3	71	0.0044	1
pir B36664 B36664	S59/4 homeotic protein - fruit ...	-3	76	0.0080	1

>sp|P27615|LIM2_RAT LYSOSOME MEMBRANE PROTEIN II (LIMP II) (85 KD LYSOSOMAL
 MEMBRANE SIALOGLYCOPROTEIN) (LGP85). >pir|A41180|A41180 74k
 lysosomal membrane protein LIMP - rat | 0.0 0.0 0.0 0.0 0.0
 >pir|JH0241|JH0241 85K lysosomal membrane sialoglycoprotein - rat |
 0.0 0.0 0.0 0.0 0.0 >gp|D10587|RATLGP85_1 LGP85 [Rattus sp.]
 >gp|M68965|RATLIMP2_1 lysosomal membrane protein [Rattus
 norvegicus]
 Length = 478

Plus Strand HSPs:

Score = 114 (55.2 bits), Expect = 1.1e-08, P = 1.1e-08
Identities = 22/64 (34%), Positives = 36/64 (56%), Frame = +2

Query: 254 LLC AVLGVVMILVMP SLIKQQV LKNVRIDPSSLSFAMWKEIPVPFYLSVYFFE VVNPS EI 433
IL + +++ V + Q + KN+ + + F W++ P+P Y+ YFF V NP EI
Sbjct: 16 LLVTSVTL LVARVFQKAVDQTI EKNMVLQNGTKVFD SWKPPLPVYIQFYFFNVTNPEEI 75

Query: 434 LRGE 445
L+GE
Sbjct: 76 LRGE 79

>pir|J01523|J01523 lysosomal membrane 85K sialoglycoprotein precursor - human |
0.0 0.0 0.0 0.0 0.0 >gp|D12676|HUMHLGP85_1 85kDa human lysosomal
sialoglycoprotein [Homo sapiens]
Length = 478

Plus Strand HSPs:

Score = 109 (52.8 bits), Expect = 6.3e-08, P = 6.3e-08
Identities = 21/64 (32%), Positives = 35/64 (54%), Frame = +2

Query: 254 LLC AVLGVVMILVMP SLIKQQV LKNVRIDPSSLSFAMWKEIPVPFYLSVYFFE VVNPS EI 433
IL + +++ V + Q + K + + + F W++ P+P Y YFF V NP EI
Sbjct: 16 LLVTSVTL LVARVFQKAVDQSI EKKIVLNGTEAFDSWKPPLPVYIQFYFFNVTNPEEI 75

Query: 434 LRGE 445
L+GE
Sbjct: 76 LRGE 79

>sp|P10284|HM26_MOUSE HOMEBOX PROTEIN HOX-2.6. >pir|A31757|A31757 homeotic
protein Hox 2.6 - mouse | 0.0 0.0 0.0 0.0 0.0 >gp|M36654|MUSHOX26_1
Mouse homeo box 2.6 (Hox-2.6) mRNA, complete cds. [Mus musculus]
Length = 250

Minus Strand HSPs:

Score = 61 (29.7 bits), Expect = 0.72, P = 0.52
Identities = 13/41 (31%), Positives = 19/41 (46%), Frame = -2

Query: 251 PRRFAPPPPSALAVPPMSRVRPGLAGLRDSRGTEATEAT 129
P P PPPP + P + V+P L G +EA ++
Sbjct: 75 PPPPPPPPPPPGLSPRAFVQPTAGALLPEPGORSEAVSSS 115

Score = 60 (29.3 bits), Expect = 2.4e-06, Poisson P(2) = 2.4e-06
Identities = 13/25 (52%), Positives = 13/25 (52%), Frame = -2

Query: 278 PHRAQRTAAPRRFAPPPPSALAVPP 204
P QR AA R P PPPP PP
Sbjct: 59 PCTVQRYAACRDPGPPPPPPPPPP 83

>sp|P16671|CD36_HUMAN PLATELET GLYCOPROTEIN IV (GPIV) (GPIIIB) (CD36 ANTIGEN).
>pir|A30989|A30989 CD36 protein - human | 0.0 0.0 0.0 0.0 0.0
>gp|M24795|HUMANTCD36_1 Human CD36 antigen mRNA, complete cds.

[Homo sapiens] >gp|M98398|HUMCD3613_1 antigen CD36 [Homo sapiens]
>gp|M98399|HUMCD3621_1 antigen CD36 [Homo sapiens]
Length = 472

Plus Strand HSPs:

Score = 94 (45.5 bits), Expect = 1.1e-05, P = 1.1e-05
Identities = 18/64 (28%), Positives = 36/64 (56%), Frame = +2

Query: 245 VVGLLCAVLGVVMILVMPSLIKQQVLKNVRIDPSSLSFAMWKEIPVPFYLSVYFFEVDNP 424
V+G + AV G +++ V LI++ + K V ++ ++F W + Y + F+V NP
Sbjct: 14 VIGAVLAVFGGILMPVGDLLIQKTIKKQVLEEGTIAFKNWVKTGTEVYRQFWIFDVQNP 73
Query: 425 SEIL 436
E++
Sbjct: 74 QEVM 77

>gp|L06850|HUMCD36B_1 antigen CD36 [Homo sapiens]
Length = 472

Plus Strand HSPs:

Score = 94 (45.5 bits), Expect = 1.1e-05, P = 1.1e-05
Identities = 18/64 (28%), Positives = 36/64 (56%), Frame = +2

Query: 245 VVGLLCAVLGVVMILVMPSLIKQQVLKNVRIDPSSLSFAMWKEIPVPFYLSVYFFEVDNP 424
V+G + AV G +++ V LI++ + K V ++ ++F W + Y + F+V NP
Sbjct: 14 VIGAVLAVFGGILMPVGDLLIQKTIKKQVLEEGTIAFKNWVKTGTEVYRQFWIFDVQNP 73
Query: 425 SEIL 436
E++
Sbjct: 74 QEVM 77

>gp|L19658|RATFAT_1 FAT gene product [Rattus norvegicus]
Length = 472

Plus Strand HSPs:

Score = 92 (44.5 bits), Expect = 2.3e-05, P = 2.3e-05
Identities = 18/65 (27%), Positives = 36/65 (55%), Frame = +2

Query: 245 VVGLLCAVLGVVMILVMPSLIKQQVLKNVRIDPSSLSFAMWKEIPVPFYLSVYFFEVDNP 424
V+G + AV G +++ V LI++ + + V ++ ++F W + Y + F+V NP
Sbjct: 14 VIGAVLAVFGGILMPVGDLLIEKTIKREVVLLEGTIAFKNWVKTGTTVYRQFWVFDVQNP 73
Query: 425 SEILK 439
E+ K
Sbjct: 74 EEVAK 78

>pir|A43932|A43932 mucin - human (fragment) | 0.0 0.0 0.0 0.0 0.0
>gp|M74027|HUMMUC2A_1 mucin [Homo sapiens]
Length = 573

Minus Strand HSPs:

JUN. 7. 2004 9:30PM

PABST PATENT GROUP

NO. 0372 P. 37

Score = 60 (29.1 bits), Expect = 1.4, P = 0.74
Identities = 12/21 (57%), Positives = 14/21 (66%), Frame = -1

Query: 279 TTPSTAHSSTTPSPTATQRP 217
TPS ++ TTPSPT T P
Sbjct: 377 TTPSPPTTMTTPSPTTTPSP 397

Score = 58 (28.1 bits), Expect = 3.8e-05, Poisson P(2) = 3.8e-05
Identities = 12/20 (60%), Positives = 14/20 (70%), Frame = -1

Query: 285 IITPSTAHSSTTPSPTAT 226
I TPS ++ TTPSPT T
Sbjct: 343 IITPSPPTTMTTPSPTT 362

>pir|B60492|B60492 homeotic protein Hox B4 - human | 0.0 0.0 0.0 0.0 0.0
Length = 251

Minus Strand HSPs:

Score = 57 (27.8 bits), Expect = 2.9, P = 0.95
Identities = 12/21 (57%), Positives = 12/21 (57%), Frame = -2

Query: 266 QRTAAPRRPAPPPPSALAVPP 204
QR AA R P PPPP PP
Sbjct: 63 QRYAACRDPGPPPPPPPPPP 83

Score = 56 (27.3 bits), Expect = 4.0e-05, Poisson P(2) = 4.0e-05
Identities = 11/20 (55%), Positives = 12/20 (60%), Frame = -2

Query: 254 APRRPAPPPPSALAVPPMSR 195
+PR PAPPP AL P R
Sbjct: 90 SPRAPAPPPAGALLPEPGQR 109

>sp|Q01200|PRIA LENEED PRIA PROTEIN. >pir|S23106|S23106 pria protein - Shiitake
mushroom | 0.0 0.0 0.0 0.0 0.0 >gp|X60956|LEPRIA_1 pria gene
product [Lentinus edodes]
Length = 258

Minus Strand HSPs:

Score = 62 (30.0 bits), Expect = 0.61, P = 0.46
Identities = 13/31 (41%), Positives = 18/31 (58%), Frame = -1

Query: 318 TCCLMSEGITRIITPSTAHSSTTPSPTAT 226
TCCL + TPS+AH + T SP++T
Sbjct: 90 TCCLPKWPTSTPTPTPSSAHHTSTHTSPSST 120

Score = 56 (27.1 bits), Expect = 5.6e-05, Poisson P(2) = 5.6e-05
Identities = 13/33 (39%), Positives = 16/33 (48%), Frame = -1

Query: 276 TPSTAHSSTTPSPTATQRPGRAAHVACAARAR 178
TPS+ +TP P+AT G H A AR
Sbjct: 143 TPSSPSKPSSTPKPSATPNKGNHGHYKRAHVAR 175

>pir|S12968|S12968 Acrosin, sperm - Pig #EC-number 3.4.21.10 | 0.0 0.0 0.0 0.0

0.0
Length = 374

Minus Strand HSPs:

Score = 59 (28.8 bits), Expect = 1.6, P = 0.79
Identities = 14/48 (29%), Positives = 24/48 (50%), Frame = -2

Query: 251 PRRPAPPPPSALAVPPMSRVRLGLRLDSRGTEATEATEATEAEAE 108
P++ + PP AL+ + ++ L G S G + TE T++ E A
Sbjct: 326 PQQVSAPPPQALSFARLQQLIEALKGTAFSSGRSYETETTDIQELPA 373

Score = 56 (27.3 bits), Expect = 6.7e-05, Poisson P(2) = 6.7e-05
Identities = 11/21 (52%), Positives = 12/21 (57%), Frame = -2

Query: 266 QRTAAPRRPAPPPPSALAVPP 204
Q + PR PAPPPP PP
Sbjct: 294 QPGSRPRPPAPPPPPPPPPPP 314

>gp|L23108|MUSCDANTI_1 CD36 antigen [Mus musculus]
Length = 473

Plus Strand HSPs:

Score = 88 (42.6 bits), Expect = 9.0e-05, P = 9.0e-05
Identities = 17/65 (26%), Positives = 35/65 (53%), Frame = +2

Query: 245 VVGLLCAVLGVVMILVMPSLIKQQVLKNVRIDPSSLSFAMWKEIPVPFYLSVYFFEVVNP 424
V+G + AV G +++ V LI++ + + V ++ + +F W + Y + F+V NP
Sbjct: 15 VIGAVLAVFGGILMPVGDMLIEKTIKREVVLEEGTAFKNWVKGTGTTVYRQFWIFDVQNP 74

Query: 425 SEILK 439
++ K
Sbjct: 75 DDVAK 79