


Exhibit 9

Lab Notebook Pages of
Dr. Susan L. Acton

Successive Rounds of screening to clone the murine SR-BI from a 3T3-L1 Adipocyte cDNA Library



Probe: 500-600 bp BamHI fragment (5') of hSR-BI in pcDNA1

3T3-L1 Adipocyte Library

pcDNA1 in MC1061/P3

G. Balclini electroporated & plated bugs, then scraped plates & froze. This is what she gave me (3 vials)

Use vial 2

Scrape a little & thaw, dilute 1:100, then 1:50 and plate 300 μ l on a 10 cm plate

She did 20 plates of 30,000 colonies (must have been very small 'colonies')

Note that the library is non-directional & is in the ~~Bst~~XI sites in pcDNA1

X158 Screen 3T3-L1 Adipocyte Library for murine SR-BI

Purpose: Clone the murine SR-BI for further analysis

Plate Library

Pour 100mm LB Amp/Tet (15/8 µg/ml) Plates
let dry 48 hrs r.t.

Scrape some bugs from frozen tube 2. (from
G. Baldwin, Lodish lab) of 3T3-L1 library in
MC1061/P3 pCDNA1 expression
thaw scraped bugs, dilute 1:100, then 1:17
and plate 50 µl

2 µl undil bugs tube #2
198 µl LB

200 µl of 1:100
↳ 118 µl of 1:100
1.88 µl of LB

2 ml of 1:17
↳ plate 50 µl/100mm plate

Incubate plates 14 hrs (6pm-8pm) 37°C
Transfer plates to 4°C /hr

Make lifts

- ① Carefully lay down sterile Nitrocellulose
Millipore HATF prenumbered, number side down
onto bugs
- ② When wet, poke 3 holes thru membrane and
plate
- ③ Lift filter off plate with blunt end forceps

- ~~Loop~~
- ④ Lay bug side up onto fresh LB A/T plate
 - ⑤ Lay another filter (pre-numbered) numby-side down onto bugs on filter to make sandwich
 - ⑥ Poke holes thru duplicate filter to match holes in original filter
 - ⑦ Let bugs grow 3 hrs 37°C between filters
 - ⑧ Let bugs regrow on original filter plate 3 hrs. Then store wrapped in parafilm 4°C.
 - ⑨ As in Maniatis, (1.98) lyse colonies
Set up 4 trays with Whatman 3MM Paper soaked with (in order)
 - ① 10% SDS
 - ② Denaturing soln. (0.5N NaOH, 1.5M NaCl)
 - ③ Neutralising soln (1.5M NaCl, 0.5M Tris-Cl pH 7.4)
 - ④ 2X SSC
 - a) Peel nitrocellulose sandwich from gel & place on SDS-impregnated 3MM paper 8-3 min
 - b) Wipe off excess with edge of tray & transfer to denat. tray 5 min
 - c) Transfer to Neutral tray 5 min
 - d) Transfer to 2X SSC tray ~~5 min~~ but separate sandwiches first - 5 min
 - e) Transfer to dry 3MM paper colony side up
Dry at least 30 min r.t.
 - ⑩ Sandwich filters between 2 sheets of dry 3MM paper. Fix by baking 1 hr 80°C in vacuum oven
 - ⑪ Store filters under vacuum no heat

X158 cont'dMake SR-BI probe

Digest #		10X Ac BSA	H ₂ O	10X NEB Buffer	enzyme
① phasRIII (1.03 µg/µl)	10 µl	5 µl	28 µl	5 µl #2	PSI/HindIII 1 µl/1 µl
② pcDNA1 (1 µg/µl)	10 µl	5 µl	28 µl	5 µl #2	↓
③ phasRIII	10 µl	5 µl	28 µl	5 µl	BamHI ↓
④ pcDNA1	10 µl	5 µl	28 µl	5 µl	buffer ↓

Incubate 37°C 2hr 2:10-4:10 freeze digest
 Ran mini-gel; all digests looked good.
 Clean up digests ① + ③

Phenol extract (50 µl)
 Back extract w/ 50 µl TE
 Precip DNA in 100 µl
 with 10 µl 3M NaOAc
 260 µl 100% EtOH
 -80°C 20 min
 Spin 30 min 4°C
 Wash w/ 70% EtOH (cold)
 Speed vac Dry
 Resuspend in 15 µl TE
 Add 5 µl 10X blue juice
 Heat 68°C 10 min

Run on 1% Low melt (Seaplaque) TAE gel

λ HindIII 0.174
 AaeIII

Digest ①
 4 lanes

Digest ③
 4 lanes


4°C overnight
 6 V
 turned up to
 20 V at noon
 Cut out bands
 at 6 pm

IX58 cont'd

1 Hind III
OX 174 Hind III

phase BI
Pst I/Hind III

phase BI
Bam HI



Probes

C. ————— A
 B

D. —————

For screening lifts, use probe B

Label Probe

- ① Melt probe B, take 10ul & add to 10ul H₂O. Label as probe B diluted 1:1
 - ② Combine 5ul diluted probe B
5ul ^{mirrors} (from Amersham Kit 9mers)
boil 15'
spin down liquid
 - ③ Add 10ul labeling buffer (kit)
23ul H₂O
2ul enzyme
5ul α -³²P dCTP
Mix by pipetting gently
Spin
Put in 37°C warm room ~45 min
2hr
- 2pm - 4pm

[X158] cont'd

Check incorporation

Combine
20 μ l salmon sperm DNA
1 μ l label rxn
180 μ l 5% cold TCA

Spin 10' 4°C
Remove 100 μ l & count as top
Count remainder as bottom

$$\% \text{ incorp} = \frac{\text{bottom} - \text{top}}{\text{bottom} + \text{top}} \times 100\% = 22\%$$

$$\% \text{ incorp} = 22\%$$

Clean up probe:

Add 5 μ l LPA to rxn
25 μ l 7.5 M NH₄ Acetate
125 μ l 100% EtOH (cold)

Put at -80°C 20 min
Spin 20 min 4°C
Rinse w/ cold 70% EtOH
Speed-vac dry
Resuspend in 30 μ l TE

Add to prehyb sol'n containing filters
Hybridize 50°C overnight

X158 cont'd

Prehyb - 50°C

500 mM PB, 1 mM EDTA, 7% SDS, 1% BSA, 100 µg/ml
salmon sperm DNA (boiled prior to adding)

For 60 mls: 30 mls 1M PB ✓
120 µl 500 mM EDTA ✓
21 mls 20% SDS ✓
0.6g BSA ✓
0.6 ml salmon sperm DNA ✓
8.28 ml H₂O ✓

Hyb - 50°C in prehyb sol'n + probe

Washes -

Wash blot 2x fast w/ 300 mM PB r.t.
wash " 1x 10' " " r.t.
Wash " 2x 10' Sol'n A (300 mM) 54°C
Wash " " " " Sol'n B (300 mM) 54°C

Put down on film - 48 hrs r.t.

- Let filters air dry on whatman
- Arrange on old piece film covered with saran wrap
- ~~Scrape~~ Tape down filters & cover w/ saran wrap
- Put marks (from Stratagene ruler) down to orient film w/ filters


Results:

> 100 small dots which are on duplicate lifts.
Background is low. Looks good!

[X158] cont'd

Secondary Lifts

- Picked 19 positive dots,
1 negative control
from primary lifts
- Used pasture p/pets, large end to pick
plugs from plates into 0.5 ml LB + Amp/Tet
overnight 37°C
- Make dilutions



- Plate bugs
First time plated way too many so replated
~13 hrs 37°C
- Transfer plates to cold room 1 hr 4°C

Lifts

- Carefully lay down sterile HATF nitrocellulose
millipore prenumbered, number side down
onto bugs
- When wet, poke 3 holes thru membrane
and plate
- Lift filter off plate with blunt end forceps
- Lay bug side up onto fresh LB A/T plate
(LB Kan for 54-60)
- Lay another filter (labeled B) number
side down onto bugs on filter to make
a sandwich
- Poke holes thru duplicate filter to match
holes in original filter

- ⑦ let bugs grow 3 hrs 37°C between filters
 ⑧ let bugs grow on original plate 3 hrs. then store wrapped in parafilm 4°C
 ⑨ As in Maniatis (1.98) lyse colonies
 Set up 4 trays with Whatman 3MM paper soaked with

3 min
 5 min
 5 min
 5 min

- ① 10% SDS
 ② Denaturing soln (0.5N NaOH, 1.5M NaCl)
 ③ Neutralizing soln (1.5M NaCl, 0.5M Tris-Cl, pH 7.4)
 ④ 2X SSC

- a) Peel nitro. sandwich from gel & place on SDS-immug. 3MM paper
 b) Wipe off excess with edge of tray & transfer to denat. tray. 5 min
 c) Transfer to neutral tray. 5 min
 d) Transfer to 2X SSC tray after separating filters
 e) Transfer to dry 3MM paper colony side up. Dry at least 30 min r.t.

- ⑩ Sandwich filters between sheets of dry 3MM paper. Fix by baking 1 hr 80°C in vacuum oven
 ⑪ Store filters under vacuum no heat

Label probe (5' BamHI fragment of hSR-BI) as for 1% HTS

ALERT 3 ID:DERENKOV PRESET TIME: 1.00 TUI 4.0
 SAMPLE REPEAT: 1 CYCLE REPEAT: 1000IN RS232C
 IN: 1.000IN RS232C
 CHANNEL 1-LL1000 UL1000 ISIGMA: 2.00 BKG STD: 0.00 BKG 2510: 0.00 LOR:
 CHANNEL 2-LL1600 UL1000 ISIGMA: 2.00 BKG STD: 0.00 BKG 7010: 0.00 LOR:
 DATA CALC: CFM, UNKNOWN REPLICATES: 1 NORM FACTOR: 01.00000
 HALF LIFE(DAYS) IN

SAN	CFM1	CFM2	TIME
1	102136.59	757.00	1.00
2	107760.97	324.00	1.00

36% incorporation

Do L...L... 11.1 - 1.0m - ~~Not a good result~~ 6 ok

Wash blots - Note: May have accidentally used 600mM PB
final rather than 300mM

- ① Wash blot 1X fast w/ 300mM PB r.t.
- ② Wash blot 1X 10' r.t. 300mM PB
- ③ Wash 2X 54°C sol'n A 300mM 10' each

Wash A: 300mM PB, 5% SDS, 0.5% BSA, 1mM EDTA

For 500ml: 300ml 0.5M PB
25g SDS
2.5g BSA
1ml 500mM EDTA
H₂O to 500ml

- ④ Wash 2X 54°C sol'n B (300mM PB) 10' each

Wash B: 300mM PB, 1% SDS, 1mM EDTA

300ml 0.5M PB
15ml 20% SDS
1ml 500mM EDTA
H₂O to 500ml

Dry & expose to film

- ① Let filters dry on whatman paper
- ② Use old piece of film as support
wrap it in saran wrap
put filters on it & tape down
put stragene ruler pieces on it to
line up films w/ filters
cover with saran wrap
- ③ Expose to film 2 days

Results: 2° lifts don't look as good as 1° - should
not have amplified plus pick from 1°. I
think high background may be due
to accidental use of 600mM PB rather
than 300mM PB.

Still there seem to be some definite
positives so picked them for 3°.

X158 cont'd

3° Lifts

- A. Picked 9 positive dots (#61-69)
1 negative control (#70)
from secondary lifts
- B. Used sterile pasture pipets, large end to
pick plugs into 0.5 ml LB. Vortexed well.

LB
plug + 500 µl

10 µl + 990 µl → plate 50 µl

10 µl + 990 µl → plate 50 µl

100 µl + 900 µl → plate 50 µl

results
way too many
way too many
just about right
~ 200/plate

Grow ~14 hrs 37°C

- C. Transferred ~~lowest~~ ^{highest} dilution plate to 4°C for lifts
- D. Performed lifts as before

Made mistake of re-incubating master plates
for ~48 hrs rather than 6-8 hrs. Plates
dried somewhat and colonies grew big

- E. Label probe as before

ALARY 6 ID:GERMKNKV PRESET TIME: 1.00 TUF
SAMPLE REPEAT: 1 CYCLE REPEAT: 1 SCRN 200.00
100.00 100.00 100.00 100.00 100.00 100.00
100.00 100.00 100.00 100.00 100.00 100.00
100.00 100.00 100.00 100.00 100.00 100.00
100.00 100.00 100.00 100.00 100.00 100.00
100.00 100.00 100.00 100.00 100.00 100.00

TIME	CPM1	CPM2
1	187661.97	212.00
2	175257.95	244.00

37.5% incorporated

Prehyb } as before 50°C, 300 mM PB
Hyb }

Wash 3° as before - 54°C, 300mM
Put down on film - w/screen

X160

Plasmid midiprep

preps: 6 single isolates from 3° screen #71-76

Day 1

1. Pick a single colony into 25 mls LB Amp/Tet and grow overnight shaking 37°C.

Day 2

2. Take 400 μ l out into freezer vial, add 100 μ l glycerol and freeze at -150°C.
3. Transfer remainder to Falcon 2059 15 ml tube on ice.
4. Spin rest in SS-34, 9000 rpm, 2 min 4°C. Add remaining half.
5. Dry pellet as much as possible.
6. Resuspend pellet in 500 μ l ice cold solution I by vigorous vortexing.
7. Add 1 ml fresh solution II (0.2 N NaOH, 1% SDS)
for 100 ml: 1 ml 2 N NaOH
0.5 ml 20% SDS
8.5 mls ddH₂O
Swirl gently until clear. Do not vortex. Leave on ice 10 min.
8. Add 750 μ l solution III (ice-cold). Close tube and mix contents by shaking vigorously several times. Store on ice 5'. A flocculent white precipitate should form.
9. Centrifuge 15', 4°C, 9000 rpm.
10. Recover supe and add equal volume of phenol:chloroform. Mix by vortexing.
11. Spin 9000 rpm, 5'.
12. Add 2 volumes, ethanol r.t., vortex, let stand 5'. Spin in SS-34 rotor for 15' 9,000 rpm.
13. Carefully aspirate off all liquid (want it dry) and rinse once in cold 70% EtOH. Aspirate again getting any extraneous drops. Airdry 10' minutes.
14. Redissolve in 100 μ l TE plus DNase-free RNase (20 μ g/ml). Vortex briefly. Incubate 37°C, 0.5 - 2 hr. Transfer to sterile eppendorf. Refrig 4°C o.n.
15. Quantitate by dilutions onto EtBr plate.
Add DNase-free RNase to final 10 μ g/ml




X160
cont'd

Digests of Midi-preps

Digest #	DNA	10X ACPSA ul	NEB Buffer 10X	H ₂ O	1ul each enzymes	DNA Prep
1	#71 - 2ul	2ul	2.2ul	14ul	Hind III	O.B.
2			2	13	Hind III / Not I	
3			2	12	Pst I	
4			2	13	Bam HI	
5			u	13		
6	#72					
7						
8						
9						
10						
11	#73					
12						
13						
14						
15						
16	#74					
17						
18						
19						
20						
21	#75					
22						
23						
24						
25						
26	#76					
27						
28						
29						
30						

37°C 11:45am - 2:45

X160 cont'd



Note that orientation B was found to be correct after sequencing of #74 + #76

X161

DEAE dextran transfections of COS M6 cells

materials:

1. 35 mm dishes
2. DMEM with 10% FBS
3. Chloroquine (40 mM in CMF PBS, sterile filtered)
4. DNA
5. CMF PBS
6. DEAE-dextran (10 mg/ml in CMF PBS (autoclaved))
7. DMSO
8. cPBS
9. sterile tips

method:

day 0 (set up cells)

Set COS M6 cells in 100 mm dishes

(Set 1 confluent T75 into 2.25 100 mm dishes; or 1.5×10^6 cells/dish in 10 ml DMEM with 10% FBS)

day 1 (transfect)

1. In sterile polypropylene tubes prepare for each dish add (in order):
 - a) add CMF PBS to 1.9 ml
 - b) DNA - $8 \mu\text{g}$ /dish
 - c) 100 μl of 10 mg/ml DEAE-dextran, vortex well (7 seconds)

tube #	plates	DNA	CMF PBS	10 mg/ml DEAE-dextran
1	1	phsR BI Cl. 1.03 mg/ml 9.7 μl	1.9 ml	100 μl
2	1	#74 (7/22/94) 1.94 mg/ml 5.2 ml	1.9 ml	100 μl
3	1	#76 (7/22/94) 0.85 mg/ml 11.8 μl	1.9 ml	100 μl
4				

2. Rinse cells with 10 mls CMF PBS (37C). Aspirate PBS.
3. Add transfection cocktail from step 1 and distribute evenly. Incubate 30 minutes, distributing liquid every 10 min.
4. Add 8 ml DMEM 10% FBS + 80 μM chloroquine and incubate 37C 2.5 hrs.
5. Aspirate off medium and replace with 5 ml 10% DMSO in DMEM 10% FBS for 2.5 min.
6. Aspirate off and wash once with 10 mls cPBS.
7. Refeed with 10 ml warm DMEM 10% FBS/dish. Incubate overnight.

X/lel cont'd

Lipoprotein 4°C Binding Assay

materials:

1. 6-well dishes
2. DMEM with 10% FBS
3. iodinated and unlabeled lipoprotein
4. Tris-HCl wash buffer
5. Tris-HCl BSA wash buffer
6. 0.1 N NaOH
7. 10 x 75 tubes

method:

day 0: Set up cells

Set transfected COS in 6-well dishes at 1×10^6 cells/well in 3 ml/well DMEM + 10% FBS + 1 mM Nabutyrate.

day 1: Binding assay

1. Cool cells down on ice for 30 minutes.
2. Refeed cells labeled ligand 1 ml/well (Hams + HEPES + 10% FBS).
 2 Hot - ^{125}I -lipoprotein: 10.4 μM ^{125}I -ALDL (#60 104 $\mu\text{g}/\text{ml}$) LDL #161 B (39 $\mu\text{g}/\text{ml}$)
 1 Hot + cold: + 46.4 μM M-BSA #230 DR (3.2 $\mu\text{g}/\text{ml}$)
 Incubate 2 hrs 4°C (in cold room) on shaker.
3. Wash cells (1 ml each)
 3X fast with BSA wash buffer
 2X 5 min with BSA wash buffer
 2X fast with Tris wash buffer
4. Add 1.0 ml 0.1 N NaOH. Leave 20 min r.t. on shaker
5. Remove 50 μl to 10 x 75 tube and freeze for Lowry.
6. Count 500 μl of rest of sample.
7. Count 10 μl of medium + label for specific activity.

#74

	-	+ LDL	+ M-BSA
1	3	5	
2	4		

#76

6	8	10
7	9	

InsR-BI

11	12	13
		14

For 14.5 ml med: 14.5 ml Hams + HEPES + 5% ALPDS
 + 145 μl ^{125}I -ALDL (#61)


5 ml
 + 51.3 μl LDL #161 B

3 ml
 + 375 μl M-BSA #230

X163 Prepare blunt-cut pcDNA1 for reversing orientation of clone #76 mSR-BI

Purpose: Reverse orientation of clone #76 which appears to be full-length mSR-BI

Approach: Cut mSR-BI out of clone #76 with Hind III/Not I ~~open~~ chew ends to blunt it, ligate into blunt & CIPed pcDNA1



Cut pcDNA1 with EcoRV which is a blunt cutter

DNA	H ₂ O	10X ACBSA	10X NEB Buff 2	EcoRV 20,000 u/ml
#76 (10.8 µg) 15.5 µl (= 20 µg) pcDNA1 (1.29 µg/µl) 10/29/92	14.5	4 µl	4 µl	2 µl

Incubate 37°C

12:50 pm - 12:50 pm (Thurs)



← after 3 hrs

Add 60 µl dd H₂O, extract w/ 100 µl phenol:chloroform:isoamylalcohol 25:24:1
Add 10 µl 3M NaOH pH=5, add 250 µl 100% EtOH, Freeze -20°C 1 hr
Spin 10 min 4°C

[X163] cont'd

4 blunt-cut pc-DNA (cont'd) - CIP

redissolve in 90 μ l 10 mM Tris-Cl (pH 8)
remove 2 μ l & save as unCIPed

To remainder add 10 μ l 10X CIP dephosph. buffer
and 8 Units of CIP:

for blunt 1 μ l / 2 pmol, estimate 20 μ g = 16 pmol
so 8 Units CIP

Incubate 15 min at 37°C 3:35 - 3:50

Then ~~add 8 Units more~~ & incubate 45 min at 55°C.
3:50 - 4:35

Kill the CIP

Add SDS & EDTA (pH 8.0) to final 0.5% & 5 mM, respect.
Mix well. \rightarrow 5 μ l 10% SDS
1 μ l 0.5 M EDTA
Add proteinase K to final 100 μ g/ml. (10 μ l of 10 mg/ml)
Incubate 30 min 56°C. (4:50 - 5:20 pm)

Cool the rxn to r.t.

Extract once w/ phenol once w/ phenol:chloroform

Add 1/10 vol 3 M NaOAc pH 5.0 (10 μ l)

Mix well add 25% EtOH (250 μ l)

Mix & store @ -20°C 2 1/2 hrs

Centrifuge

wash pellet w/ 70% EtOH

Redissolve in TE

Ran overnight on 1% Nusieve (should have been Seaplaque)
Agarose gel. Band was wavy, but cut it out anyway

X164 Prepare #76 mSR-BI cDNA for reinserting into pCDNA1 in correct orientation


Purpose: Cut out #76 insert to be reinserted in sense orientation. #76 plasmid appears to be full-length mSR-BI but is in the incorrect orientation in pCDNA1.

Digest #76 plasmid:

	DNA	H ₂ O	10X AcBSA	10X NotI buffer	HindIII	NotI
①	#76 (0.85 µg/µl) 24 µl (= 20 µg)	4 µl	10	10 (NotI)	5	5
②	#76 24 µl	4 µl	10	10 µl (NEB L)	5	5

5:45 pm 8:15 pm 37°C

Remove 2 µl from each and analyze on 1% SeaKem GT6 Agarose TAE gel



To the remainder add 10 µl 3M NaAc & 250 µl 100% EtOH
Store at -20°C for 15 min
Spin 10 min 4°C
Wash w/ 70% EtOH
Resuspend in TE (30 µl)
Add 4 µl 10X blue juice

To blunt ends, add 5 μ l 10X T4 DNA Polymerase Buffer
 + 2 μ l of dNTP mix (25mM) for final 1mM
 + 7 μ l H_2O = 50 μ l rxn ~~37° 30 min~~
 + 1 μ l T4 DNA Polymerase 3U/ μ l 37° 30 min (11:30am-12 noon)

Add 4 μ l 0.25M EDTA
 Add 55 μ l Phenol/Chloroform/Isoamyl Alcohol (25:24:1)
 Shake, spin, save aq. layer

To aq. phase
 Add 10

1/10 V 5M NaCl (= 5 μ l)
 Shake add 5 μ l of 1X linear polyacryl.
 add 150 μ l 100% EtOH
 Chill @ -80°C 10'
 Spin 10' 4°C
 Wash pellet w/ 70% EtOH
 Speed vac dry

Resuspend pellet in 36 μ l TE
 Run 2 μ l on mini-gel



Add blue juice & run on 1% TAE SeaPlaque agarose gel
 in cold room at 20 V (med size app)
 5:45pm - 11:00 am
 Cut out insert (~2.3 kb) & ligate into vector (see following page.)

1X16.5

Ligation Rxns

- 1) Run blunt vector + blunt insert on 1% GTC Seakem gel to quantitate relative amounts

vector
DNA insert



Estimate that 3 μ l vector = 3 μ l cDNA

Want ~ twice as much insert as vector
in 10 μ l rxn so use
4 μ l vector and 6 μ l insert

- 2) Heat tubes to 70°C for 10-15 min to melt agarose
3) Combine aliquots of melted gel slices in tubes prewarmed to 70°C

tube	vector		blunt cDNA insert	H ₂ O	2X ligation rxn mix (10-100 μ l)
	CTP	not CTP			
1	4 μ l	-	-	6 μ l	10 μ l
2	4 μ l	-	6 μ l	4 μ l	↓
3	4 μ l	-	6 μ l	4 μ l	↓
4	-	4 μ l	-	6 μ l	↓

2X Lig rxn mix:

10 μ l	1M Tris, Cl ⁻	pH 7.5
1 μ l	1M MgCl ₂	
1 μ l	1M DTT	
1.2 μ l	80 mM ATP	
86.3 μ l	H ₂ O	
100 μ l		

Take 38 μ l & add 4 μ l T4 DNA Ligase

Incubate 48 hr 12-16°C (2 pm Sat - 17 pm Mon)

X/65 cont'd

Transformation of Ligations into MC1061/P3 E. coli

1. Melt the agarose ligation mixtures at 70°C for 10-15 min.
2. Meanwhile thaw competent MC1061/P3 (from J. Ashkenazi)

transform				Results # colonies
1	Lig #1	CIPed vector	5µl	56
2	Lig #2	insert (cDNA) only	5µl	2.8
3	Lig #3	CIPed vector + insert	5µl	30
4	Lig #4	uncIPed vector	5µl	>500
5	TE only		5µl	47
6	pcDNA4 - 100ng/µl		1µl	>500

3. Add 5µl of melted agarose ligation mixtures or control aliquot.
Mix by gentle shaking
Store on ice 30 min 12:45-1:15
4. Transfer to circulating water bath at 42°C
Incubate exactly 90 sec. DO NOT SHAKE TUBES.
5. Rapidly transfer to ice bath. Allow cells to chill 1-2 min.
6. Add 400µl LB medium warm to 37°C & transfer to shaking incubator at 37°C. Incubate 1 hr
1:25-2:25
7. Plate onto LB Amp/Tet (15µg/µl/8µg/µl) plates
100µl/100mm plate
Incubate 37°C overnight

Results: Insert did not ligate into CIPed vector. Vector is CIPed well. Try ~~the~~ blunting cDNA again and GeneClean both it and vector.

X116

GeneClean = vector + cDNA insert (#76)

	volume	^{2.5V} NaI	glassmilk
insert -	150 μ l	375 μ l	10 μ l
vector -	150 μ l	375 μ l	10 μ l

Since ligation did not work between vector + cDNA,
need to make sure ends of cDNA are blunt


Remove protruding 3' termini from cDNA insert

Combine: 15 μ l of genecleaned mSR-BI cDNA insert
 1 μ l of 1:1:1:1 mix of dNTPs (25mM of each) \Rightarrow final 1.25mM
 2 μ l of NEB Buffer 2 10X (final: 10mM Tris, 10mM MgCl₂, 50mM NaCl,
 1mM DTT pH 7.9)
 2.0 μ l ^{add H₂O} T4 DNA Polymerase (3 U/ μ l)
 Incubate 15 min 12°C

Inactivate polymerase by heating to 75°C for 10 min (11-11:10 am)

Run vector + cDNA insert on minigel to quantitate amts

1 μ l + 2 μ l TE $\xrightarrow{1 \mu l}$ + 2 μ l TE
 ① $\xrightarrow{1 \mu l}$ ②
 run 3 μ l + 2 μ l 2x blue juice = 0.67
 run 3 μ l + 2 μ l 2x blue juice = 0.22



For a 1 vector : 2 insert molar ratio
 use equal μ g of vector + insert
 Looks like vector stock = 1.5x cDNA

Ligate cDNA into vector

Ligation	blunt pcDNA1 CIPed	msr-BI blunt cDNA	H ₂ O	BRL 5 10X T4 DNA Lig buffer	BRL T4 DNA Ligase (4U/μl)
1	3μl	-	4.5μl	2μl	0.5μl
2	3μl	4.5μl	-	2μl	0.5μl
3	1μl unCIPed vector -	-	6.5μl	2μl	0.5μl

Mix well, spin & incubate 12-16°C (~18°C) 12:15pm - 10:30 PM
(~48 hrs)

Transform MC1061/P3

transf	amt	results # colonies
1 Lig #1	5μl	123
2 Lig #2	5μl	32
3 Lig #3	5μl	33
4 TE only	5μl	
5 pcDNA1-3000μg/ml 4μl	4μl	>1000

transformation worked but ~~some~~ but other results don't make sense

- 1) Add 5μl of above DNA to thawed competent MC1061/T3.
Mix by gentle shaking
Store on ice 30' (10:40-11:10)
- 2) Transfer to circulating water bath at 42°C.
Incubate exactly 90 sec.
- 3) Rapidly transfer to ice bath. Allow cells to chill 1-2 min.
Add 100μl LB. Warm to 37°C & transfer to shaking incubator at 37°C.
- 4) Plate onto LB Amp^r (15μg/ml / 8μg/ml)
100μl / 100mm dish Incubate 37°C O.N.

Run ligation rxns on gel to check



Lig



- #1 - one piece but 6 Kb! Should be 4 Kb. Must have cut the wrong vector? Redo preparation of blunt-cut cIPed pcDNA1
- #2 - some ligations appear to have worked possibly cDNA - cDNA worked best. Next time cut way back on amt of cDNA: vector.
- #3 - can't see any DNA

see below (1.5 Kb)

Run with 1 Kb standards

Results:

blunt cut pcDNA1 is 4 Kb as it should be
XbaI stds are junk



X167

Redo Ligations

<u>Ligation</u>	blunt pcDNA1 CIPed	MSR-B1 blunt cDNA	H ₂ O	BRL 5X T4 DNA Lig buffer	BRL T4 DNA Ligase (14 μ l)
1	1 μ l	-	14	4 μ l	1 0 μ l
2	1 μ l	1 μ l	13	4	1
3	1 μ l	1 μ l of 1:5	13	4	1
4	1 μ l of 1:5	1 μ l of 1:5	13	4	1
5	1 μ l of 1:5	-	14	4	1

Mix well, spin & incubate 17-16°C (submerged) O.N. (5pm - 2pm)

Transform

			Results # colonies
1	Lig #1 - 5 μ l	780 μ l comp MC1061/TS CRO	68
2	2 - 5 μ l		136
3	3 - 5 μ l		136
4	4 - 5 μ l		65
5	5 - 5 μ l		19 + ~30 in 5 min longer
6	TE - 5 μ l		5
7	pcDNA1 - 1 μ l		72000

On ice 30' 2:30 - 3pm
Shake 42°C 90 sec

On ice 1-2 min

Add 40 μ l LB & shake

Plate onto LB A/T (37°C 45') (#1-4 my plates, #5-7 X-ray's plates)

100 μ l

PCR colonies

Purpose: Determine if colonies from transformation plate #2 contain mR-BI in correct orientation.

Set A (determine if insert within): T7 + ST6 primers

Set B (determine orientation): T7 + 05A4.4B

For 55 rxns of 70 μ l each

110 μ l 10x PCR buffer (GG)
 275 μ l T7 primer (20 μ mol/l) (VA)
 725.7 μ l H₂O
 2.2 μ l dATP
 2.2 μ l dCTP
 2.2 μ l dGTP
 2.2 μ l dTTP
 5.5 μ l Tag polymerase

(A)	412.5 μ l	412.5 μ l (B)
	+ 137.5 μ l ST6 (20 μ mol/l)	+ 2 μ l stock 05A4.4B (unknown conc.) + 135.5 μ l H ₂ O

Pick 20 colonies


#21 neg control colony from transformation plate #1
 #22 control #76 pure plasmid in SR BI in pCDNA1 backward

From 20 random colonies

10 will probably be vector only (no insert): background...
 10 will " contain DNA
 of which 5 will be backwards
 and 5 will be correct orientation

Run gel on PCR products

Run 1% GTC-Agarose. 1X TAE



Results: Only the positive control pure plasmid #76 (pcDNA1/M⁺ mSRβ1⁺ mck/m⁺) worked with the T7 + 5P6 primers. However, a band of ~1.3 Kb was present in #10, #14, & #17. The expected size is 1.35 Kb for mSRβ1 in correct orientation. Hopefully these are real.

Plasmid midiprep

preps:

2, # 10, # 14, # 17 colonies from ligation/transform plate 2.

Day 1

1. Pick a single colony into 25 mls LB Amp/Tet and grow overnight shaking 37°C.

Day 2


2. Take 400 µl, put into freezer vial, add 100 µl glycerol and freeze at -150C.
3. Transfer half to Falcon 2059 15 ml tube on ice.
4. Spin in SS-34, 9000 rpm, 2 min 4C. Discard supe, add remaining bacteria and respin.
5. Dry pellet as much as possible.
6. Resuspend pellet in 500 µl ice cold solution I by vigorous vortexing.
7. Add 1 ml fresh solution II (0.2 N NaOH, 1% SDS)
for 100 ml: 1 ml 2 N NaOH
0.5 ml 20% SDS
8.5 mls ddH2O
Swirl gently until clear. Do not vortex. Leave on ice 10 min.
8. Add 750 µl solution III (ice-cold). Close tube and mix contents by shaking vigorously several times. Store on ice 5'. A flocculent white precipitate should form.
9. Centrifuge 15', 4°C, 9000 rpm.
10. Recover supe and add equal volume of phenol:chloroform. Mix by vortexing.
11. Spin 9000 rpm, 5'. Recover aq. phase.
12. Add 2 volumes, ethanol r.t., vortex, let stand 5'.
Spin in SS-34 rotor for 15' 9,000 rpm.
13. Carefully aspirate off all liquid (want it dry) and rinse once in cold 70% EtOH. Aspirate again getting any extraneous drops. Airdry 10' minutes.
14. Redissolve in 100 µl TE plus DNase-free RNase (20 µg/ml). Vortex briefly. Incubate 37C (0.5 - 2 hr).
Transfer to sterile eppendorf.
15. Quantitate by dilutions onto EtBr plate.

Digests

Digest	DNA	10X	10X	10X	10X
1	colony 10	2ul	2ul	2ul	2ul
2					
3					
4					
5	colony 10				
6					
7					
8					
9					
10	colony 14				
11					
12					
13	colony 17				
14					
15					
16					
17	#76 msp.B1 (backwards)				
18					
19					
20					

Incubate 57°C
12:45 - 3:45

Run on 1% Seakem GTG TAE gel



Results: Neg control #2 looks like
pcDNA1 only. #10 #14 #17
all look same with insert.
size ~ 2.5 Kb. Not identical
to #76 msp.B1 (backwards).
Everything looks good!
I think all three
#10 #14 #17 are all
in correct direction.
see next page

X168

DEAE dextran transfections of COS M6 cells**materials:**

- | | |
|---|--|
| 1. 35 mm dishes | 5. CMF PBS |
| 2. DMEM with 10% FBS | 6. DEAE-dextran (10 mg/ml in CMF PBS (autoclaved)) |
| 3. Chloroquine (40 mM in CMF PBS, sterile filtered) | 7. DMSO |
| 4. DNA | 8. cPBS |
| | 9. sterile tips |

method:**day 0 (set up cells)**

Set COS M6 cells in 100 mm dishes

(Set 1 confluent T75 into 2.25 100 mm dishes; or 1.5×10^6 cells/dish in 10 ml DMEM with 10% FBS)**day 1 (transfect)**

1. In sterile polypropylene tubes prepare for each dish add (in order):

a) add CMF PBS to 1.9 ml

b) DNA - 5 μ g/dishc) 100 μ l of 10 mg/ml DEAE-dextran, vortex well (7 seconds)

tube	# plates	DNA		CMF PBS	10 mg/ml DEAE-dextran
1	1	pcDNA1 (1.29 μ g/ μ l)	15 μ g	3.9 ml	1.9 ml
2	3	pmsp-BL (0.11 μ g/ μ l)	15 μ g	30 μ l	5.7 ml
3					100 μ l
4					300 μ l

2. Rinse cells with 10 mls CMF PBS (37C). Aspirate PBS.
3. Add transfection cocktail from step 1 and distribute evenly. Incubate 30 minutes, distributing liquid every 10 min. 1:30 - 2:00
4. Add 8 ml DMEM 10% FBS + 80 μ M chloroquine and incubate 37C 2.5 hrs. 4:30
5. Aspirate off medium and replace with 5 ml 10% DMSO in DMEM 10% FBS for 2.5 min.
6. Aspirate off and wash once with 10 mls cPBS.
7. Refeed with 10 ml warm DMEM 10% FBS/dish. Incubate overnight.

Lipoprotein 4°C Binding Assay**materials:**

1. 6-well dishes
2. DMEM with 10% FBS
3. iodinated and unlabeled lipoprotein
4. Tris-HCl wash buffer
5. Tris-HCl BSA wash buffer
6. 0.1 N NaOH
7. 10 x 75 tubes

method:**day 0: Set up cells**

Set transfected COS in 6-well dishes at 1×10^6 cells/well in 3 ml/well DMEM + 10% FBS + 1 mM Nbutyrate.

day 1: Binding assay


1. Cool cells down on ice for 30 minutes.
2. Refeed cells labeled ligand 1 ml/well (Hams + HEPES + 5% NCLPDS + 10% FBS).
2 Hot - ^{125}I -lipoprotein: 10 ng/ml ^{125}I -AcLDL #62 (0.55 $\mu\text{Ci}/\mu\text{g}$)
+ Hot + cold: competitor (at 400 ng/ml)
Incubate 2 hrs 4°C (in cold room) on shaker.
3. Wash cells (1 ml each)
3X fast with BSA wash buffer
~~2X 5 min with BSA wash buffer~~
2X fast with Tris wash buffer
4. Add 1.0 ml 0.1 N NaOH. Leave 20 min r.t. on shaker
5. Remove 50 μl to 10 x 75 tube and freeze for Lowry.
6. Count 500 μl of rest of sample.
7. Count 10 μl of medium + label for specific activity.

AcLDL #55 4.6 ng/ml
LDL #16/B 39
HDL #3 31.2
VLDL #3 26.7

309 μl ^{125}I -AcLDL #62 + 17 ml Hams + 5% NCLPDS + 10 mM HEPES

4 ml + 348 μl AcLDL #55
2 ml + 20.5 μl LDL #16/B
2 ml + 25.6 μl HDL #3
2 ml + 38.6 μl VLDL #3

msr-BI in correct orientation:



msr-BI (#76) in incorrect orientation:

